

TECHNICAL BULLETIN



19 Motivation Dve, Wangara, WA, 6065 AUSTRALIA
T +61 8 9302 4000 | FREE 1800 999 196 | F +61 8 9302 5000

QT 200 CONCENTRATED DISINFECTANT

MATERIAL & FUNCTION

QT 200 is a highly active sanitizer, deodorant and fungicide based on quaternary ammonium compounds, designed for use in areas of food preparation.

QT 200 is suitable for use in food and meat export establishments.

APPEARANCE

Clear colourless liquid with distinctive odour.

APPLICATION

QT 200 is suitable for the disinfection of various types of plant and equipment, namely homogenizers, pans, tanks, bench tops, bowls, refrigerators, heat exchangers, trays, conveyor belts and mixing, mincing and kneading plants.

For non-rinse applications where food may come in contact, quat concentrations must not exceed 200 ppm (1 mL **QT200** per L). **QT 200** is recommended for fogging, spraying or manual applications at use dilutions of 0.75 mL (150 ppm quat) to 1.5 mL (300 ppm quat) with one litre of water.

QT 200 solutions should only be applied to pre-cleaned surfaces. Residues shall be removed from the product contacting surfaces of equipment, by thorough rinsing with clean water before the commencement of operations, but need not be removed from the outside equipment or from walls, floors or ceiling unless their presence is considered hazardous by the establishment officer-in-charge of inspection.

QT 200 is safe to use on all painted and metal surfaces.

PACKAGING

5, 15 & 200 Litre containers.

IMPORTANT NOTICE TO CUSTOMER

*Since the use of this product is beyond the control of either seller or manufacturer, their only obligation shall be to replace any quantity of product which is proven defective. They cannot assume any risk or liability in excess of the purchase price of the product itself, which does not include labour or any consequential damages resulting from the use of this product. Determining the suitability of this product for any intended use shall be solely the responsibility of the user. **ALWAYS TEST FIRST.***

[F:\TECHBULL\BLANK.DOC](#)